

 **Netti**® MOBILE FAMILY  
COMFORT

**Netti**®  
BY ALU REHAB



TWO  
NEW  
MODELS



*Enable joy of life*

[My-Netti.com](http://My-Netti.com)

# AN INNOVATIVE MULTIFUNCTIONAL POWER CHAIR

THE NETTI MOBILE SETS COMPLETELY  
NEW STANDARDS FOR COMFORT IN POWERED WHEELCHAIRS:

- Netti Seating Solution for comfortable, healthy and stable seating
- Optimal solutions used on other Netti wheelchairs including ergonomical backrests and adjustable foot supports and extension pieces
- A full range of motorized seat functions
- Control systems: Dynamic Shark and DX2
- Compact and highly maneuverable
- Easy to drive



## PERSONALIZATION

Netti Mobile is highly adaptable to meet the needs of a very wide range of people including those with reduced mobility, people who need help with movement or people with cognitive difficulties.

Please don't hesitate to contact us for advice in choosing the right solution.

## UNIVERSAL DEVICE RACK

With adjustable mounting for the secure attachment of a wide range of devices including oxygen concentrators, saliva suction devices.



# INTRODUCING THE NETTI MOBILE INTUITIVE

NETTI MOBILE INTUITIVE SETS THE STANDARD FOR ATTENDANT CONTROLLED POWER CHAIRS: JUST LIKE PUSHING A MANUAL CHAIR EXCEPT WITHOUT THE EFFORT.

- Automatic speed adjustment
- Easy to operate
- No joystick to maneuver
- Ergonomic design
- Fully adjustable for comfortable driving

NEW  
MODEL



- + The intuitive attendant control has a simple dashboard with an on/off switch, a button to change between slow mode for driving in tight places, and a switch for the lights in case you want to drive outside at night.
- + Press the on/off switch, wait until the green LED stops flashing, then grip the control bar and you are ready to go.
- + Just push or pull the chair in the direction you want to go as you would with a manual chair. It will drive smoothly in the direction and at the speed you require with the minimum of effort.
- + When you want to stop, just let go of the control bar.
- + Netti Mobile Intuitive has the same powerful motors and full sized wheels as the Netti Mobile. That means it can be driven outdoors without the compromises that apply to most attendant only motorised wheelchairs.
- + All this and the class leading Netti Seating Solution that for a generation has defined what is meant by a multifunctional seating system.

# THE FIRST **DYSTONIA** **ADAPTED** POWER WHEELCHAIR

Netti Dynamic Mobile is an advanced mobility aid for disabled users affected by dystonia with extension movements resulting in very strong spasms and contractions which affect the user's muscle-bone apparatus causing joint dislocations and challenging the strength of the wheelchair.

The patented Netti Dynamic System accommodates the user's extension movements letting the wheelchair work in synergy with the user's movements. When accommodating the spastic extension pattern, muscle tone and frequency of spasms will be reduced. The wheelchair is dynamic and will follow the users movement of both upper and lower body.

- Leg movements
- Hip movements
- Back movements
- Head movements
- Trunk movements



## **Netti**® CUSHIONS

FOR A COMFORTABLE, HEALTHY  
AND STABLE SITTING POSITION




 **Netti** Kyphotic



 **Netti** Stabil



 **Netti** Super Stabil



 **Netti** Uno Back



 **Netti** Smart



 **Netti** Sit



 **Netti** Uno Seat



 **Netti** Dry



 **Netti** Contour



**VITAL BASE** Venti



**VITAL BASE** Royal



# A POWER WHEELCHAIR FOR EVERY NEED



## NETTI MOBILE

The first powered wheelchair with the class leading Netti Seating System.



## NETTI MOBILE INTUITIVE

The attendant only multifunctional powered wheelchair with the Netti intuitive controller.



## NETTI DYNAMIC MOBILE

The first power wheelchair with dynamic multifunctional seating designed specifically for users with Dystonia.

## BASIC CONFIGURATIONS

### Manually adjustable

- Seat tilt in steps (-4°, -8°, -12°).
- Continuously adjustable back rest recline with gas spring.
- Netti Mobile or Netti Dynamic Mobile - Dynamic Shark control system.
- Netti Mobile Intuitive - Intuitive attendant control.

### Electric Seat Tilt

- Continuously adjustable back rest recline with gas spring.
- Netti Mobile or Netti Dynamic Mobile - Dynamic Shark control system.
- Netti Mobile Intuitive - Intuitive attendant control.

### Electric Seat Tilt & Backrest Recline

- Netti Mobile or Netti Dynamic Mobile - Dynamic Shark control system.
- Netti Mobile Intuitive - Intuitive attendant control.

### Electric Seating Tilt, Backrest Recline & Legrests

- Netti Mobile or Netti Dynamic Mobile - Dynamic DX2 expandable control system.
- Netti Mobile Intuitive - Intuitive attendant control.

- ⊕ Easy to drive and comfortable to use
- ⊕ Rated for indoor and outdoor use
- ⊕ The Netti Seating Solution makes it ideal for users who have to sit in a chair for long periods
- ⊕ Comprehensive range of basic configurations
- ⊕ Wide range of accessories and special adaptations

FEEL FREE TO CONTACT US:



Alu Rehab AS  
Bedriftsvegen 23  
N-4353 Klepp Stasjon  
Norway

T: +47 51 78 62 20  
post@My-Netti.com

Alu Rehab ApS  
Kløftehøj 8  
DK-8680 Ry  
Denmark

T: +45 87 88 73 00  
F: +45 87 88 73 19  
info@My-Netti.com

[My-Netti.com](https://www.my-netti.com)